



ZSCA2 Polyclonal Antibody

Catalog No	YP-Ab-06519
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ZSCAN2 ZFP29 ZNF854
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZSCA2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZNF854;Zinc finger protein 854;Zinc finger protein 29 homolog;Zinc finger and SCAN domain-containing protein 2;Zfp-29;ZSCA2;ZSCAN2 ZFP29 ZNF854;Zinc finger and SCAN domain-containing protein 2 (Zinc finger protein 29 homolog) (Zfp-29) (Zinc finger protein 854)
Observed Band	70 kDa
Calculated Molecular Weight	70 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Trachea,
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:May be involved in transcriptional regulation during the post-meiotic stages of spermatogenesis.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 SCAN box domain.,similarity:Contains 14 C2H2-type zinc fingers.,
Background	The protein encoded by this gene contains several copies of zinc finger motif, which is commonly found in transcriptional regulatory proteins. Studies in mice show that this gene is expressed during embryonic development, and



specifically in the testis in adult mice, suggesting that it may play a role in regulating genes in germ cells. Alternative splicing of this gene results in several transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images