



HOMEZ Polyclonal Antibody

Catalog No	YP-Ab-06525
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	HOMEZ KIAA1443
Immunogen	Synthesized peptide derived from human protein . at AA range: 280-360
Specificity	HOMEZ Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HOMEZ;KIAA1443;HOMEZ KIAA1443;Homeobox and leucine zipper protein Homez (Homeodomain leucine zipper-containing factor)
Observed Band	65 kDa
Calculated Molecular Weight	525 aa, 59 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.
Function	function:May function as a transcriptional regulator.,similarity:Contains 2 homeobox DNA-binding domains.,subunit:Homodimer or heterodimer (Potential) . Interacts with HOXC8.,tissue specificity:Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.,
Background	function:May function as a transcriptional regulator.,similarity:Contains 2 homeobox DNA-binding domains.,subunit:Homodimer or heterodimer (Potential) . Interacts with HOXC8.,tissue specificity:Ubiquitous. Strongly expressed in adult testis and kidney as well as fetal lung and kidney.,

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images