



ASPM Polyclonal Antibody

Catalog No	YP-Ab-06533
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	ASPM MCPH5
Immunogen	Synthesized peptide derived from human protein . at AA range: 1230-1310
Specificity	ASPM Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-p 1:50-300. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ASP;Asp homolog;ASPM;Calmbp1;MCPH5;ASPM MCPH5;Abnormal spindle-like microcephaly-associated protein (Abnormal spindle protein homolog) (Asp homolog)
Observed Band	382kD
Calculated Molecular Weight	382kD
Cell Pathway	Cytoplasm . Cytoplasm, cytoskeleton, spindle . Nucleus . The nuclear-cytoplasmic distribution could be regulated by the availability of calmodulin (By similarity). Localizes to spindle poles during mitosis (PubMed:19690332). Associates with microtubule minus ends (By similarity). .
Tissue Specificity	Colon adenocarcinoma,Epithelium,Fetal brain,Kidney,Lymph,Tongue,
Function	disease:Defects in ASPM are the cause of microcephaly primary type 5 (MCPH5) [MIM:608716]; also known as true microcephaly or microcephaly vera. Microcephaly is defined as a head circumference more than 3 standard deviations below the age-related mean. Brain weight is markedly reduced and the cerebral cortex is disproportionately small. Despite this marked reduction in size, the gyral pattern is relatively well preserved, with no major abnormality in cortical architecture. Primary microcephaly is further defined by the absence of other syndromic features or significant neurological deficits. This entity is



inherited as autosomal recessive trait.,function:Probable role in mitotic spindle regulation and coordination of mitotic processes (By similarity). May have a preferential role in regulating neurogenesis.,similarity:Contains 2 CH (calponin-homology) domains.,similarity:Contains 39 IQ d

Background

This gene is the human ortholog of the *Drosophila melanogaster* 'abnormal spindle' gene (*asp*), which is essential for normal mitotic spindle function in embryonic neuroblasts. Studies in mouse also suggest a role of this gene in mitotic spindle regulation, with a preferential role in regulating neurogenesis. Mutations in this gene are associated with microcephaly primary type 5. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, May 2011],


matters needing attention

Avoid repeated freezing and thawing!


Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Immunohistochemical analysis of paraffin-embedded human spleen tissue. 1,primary Antibody was diluted at 1:200(4° overnight). 2, Sodium citrate pH 6.0 was used for antigen retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200



Western Blot analysis of K-562 cell ,using primary antibody at 1:1000 dilution. Secondary antibody(catalog#:RS23920) was diluted at 1:10000