



EBP2 Polyclonal Antibody

Catalog No	YP-Ab-06558
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	EBNA1BP2 EBP2
Immunogen	Synthesized peptide derived from human protein . at AA range: 140-220
Specificity	EBP2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EBNA1 binding protein 2;EBNA1BP2;EBP2;NOBP;Nucleolar protein p40;P40;EBNA1BP2 EBP2;Probable rRNA-processing protein EBP2 (EBNA1-binding protein 2) (Nucleolar protein p40)
Observed Band	33kD
Calculated Molecular Weight	33kD
Cell Pathway	Nucleus, nucleolus . Associated with the nucleolus in an RNA-dependent manner.
Tissue Specificity	Ubiquitous.
Function	function:Required for the processing of the 27S pre-rRNA.,similarity:Belongs to the EBP2 family.,subunit:Specifically interacts with EBV EBNA1. The EBNA1-EBP2 interaction is important for the stable segregation of EBV episomes during cell division.,tissue specificity:Ubiquitous.,
Background	function:Required for the processing of the 27S pre-rRNA.,similarity:Belongs to the EBP2 family.,subunit:Specifically interacts with EBV EBNA1. The EBNA1-EBP2 interaction is important for the stable segregation of EBV episomes during cell division.,tissue specificity:Ubiquitous.,


**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night