



ZNF71 Polyclonal Antibody

Catalog No	YP-Ab-06578
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ZNF71 EZFIT
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZNF71 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EZFIT;zinc finger protein 71;ZNF71;ZNF71 EZFIT;Endothelial zinc finger protein induced by tumor necrosis factor alpha (Zinc finger protein 71)
Observed Band	54 kDa
Calculated Molecular Weight	489 aa, 55 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.
Function	function:May be involved in transcriptional regulation.,induction:By tumor necrosis factor alpha (TNFalpha).,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.,
Background	function:May be involved in transcriptional regulation.,induction:By tumor necrosis factor alpha (TNFalpha).,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 13 C2H2-type zinc fingers.,tissue



specificity:Highly expressed in placenta, followed by brain, testis, pancreas, heart, small intestine, muscle, uterus, prostate and peripheral blood leukocytes. Not detected in liver, lung, colon, stomach, salivary and thyroid gland.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images