



TF2L1 Polyclonal Antibody

Catalog No	YP-Ab-06581
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TFCP2L1 CRTR1 LBP9
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TF2L1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRTR1;CP2-related transcriptional repressor 1;Transcription factor LBP 9; Transcription factor CP2-like protein 1;LBP9;TF2L1;TFCP2L1 CRTR1 LBP9; Transcription factor CP2-like protein 1 (CP2-related transcriptional repressor 1) (CRTR-1) (Transcription factor LBP-9)
Observed Band	55 kDa
Calculated Molecular Weight	55 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in placental JEG-3 cells and very low levels of expression in non-steroidogenic cells. No expression was seen in adrenal NCI-H295A cells or in adrenal tissue.
Function	function:Transcriptional suppressor. May suppress UBP1-mediated transcriptional activation. Modulates the placental expression of CYP11A1., similarity:Belongs to the grh/CP2 family. CP2 subfamily.,tissue specificity:Highly expressed in placental JEG-3 cells and very low levels of expression in non-steroidogenic cells. No expression was seen in adrenal NCI-H295A cells or in adrenal tissue.,

**Background**

function:Transcriptional suppressor. May suppress UBP1-mediated transcriptional activation. Modulates the placental expression of CYP11A1., similarity:Belongs to the grh/CP2 family. CP2 subfamily.,tissue specificity:Highly expressed in placental JEG-3 cells and very low levels of expression in non-steroidogenic cells. No expression was seen in adrenal NCI-H295A cells or in adrenal tissue.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images