




FAIM1 Polyclonal Antibody

Catalog No	YP-Ab-06607
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	FAIM FAIM1
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	FAIM1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAIM1;FAIM FAIM1;Fas apoptotic inhibitory molecule 1
Observed Band	19kD
Calculated Molecular Weight	19kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain,Testis,
Function	function:Plays a role as an inducible effector molecule that mediates Fas resistance produced by surface Ig engagement in B cells.,similarity:Belongs to the FAIM1 family.,
Background	Fas apoptotic inhibitory molecule(FAIM) Homo sapiens The protein encoded by this gene protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night