



HS12A Polyclonal Antibody

Catalog No	YP-Ab-06615
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	HSPA12A KIAA0417
Immunogen	Synthesized peptide derived from human protein . at AA range: 230-310
Specificity	HS12A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Heat shock 70 kDa protein 12A;Heat shock protein family A member 12A; KIAA0417;HS12A;HSPA12A KIAA0417
Observed Band	75 kDa
Calculated Molecular Weight	675 aa, 75 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Widely expressed with highest levels in brain, kidney and muscle.
Function	similarity:Belongs to the heat shock protein 70 family.,tissue specificity:Widely expressed with highest levels in brain, kidney and muscle.,
Background	similarity:Belongs to the heat shock protein 70 family.,tissue specificity:Widely expressed with highest levels in brain, kidney and muscle.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

