




I20L2 Polyclonal Antibody

Catalog No	YP-Ab-06623
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ISG20L2 HSD-38 HSD38
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	I20L2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.1.-.-;HSD-38;Interferon-stimulated 20 kDa exonuclease-like 2;I20L2; ISG20L2 HSD-38 HSD38;Interferon-stimulated 20 kDa exonuclease-like 2 (EC 3.1.-.-)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	Brain,Testis,
Function	function:3'-> 5'-exoribonuclease involved in ribosome biogenesis in the processing of the 12S pre-rRNA. Displays a strong specificity for a 3'-end containing a free hydroxyl group.,similarity:Contains 1 exonuclease domain.,
Background	This gene encodes a 3'-5' exoribonuclease that may be involved in the processing of the 12S pre-rRNA. Pseudogenes have been identified on chromosomes 6 and 11. [provided by RefSeq, Dec 2014],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night