




IGFL1 Polyclonal Antibody

Catalog No	YP-Ab-06629
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IGFL1 UNQ644/PRO1274
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50
Specificity	IGFL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IGF like family member 1;IGFL1;UNQ644;IGFL1 UNQ644/PRO1274;Insulin growth factor-like family member 1
Observed Band	12kD
Calculated Molecular Weight	12kD
Cell Pathway	Secreted .
Tissue Specificity	Detected in ovary and spinal cord.
Function	similarity:Belongs to the IGFL family.,
Background	IGF like family member 1(IGFL1) Homo sapiens The protein encoded by this gene is a member of the insulin-like growth factor family of signaling molecules. The encoded protein is synthesized as a precursor protein and is proteolytically cleaved to form a secreted mature peptide. The mature peptide binds to a receptor, which in mouse was found on the cell surface of T cells. Increased expression of this gene may be linked to psoriasis. [provided by RefSeq, Aug 2016],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night