



KT33A Polyclonal Antibody

Catalog No	YP-Ab-06644
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	KRT33A HHA3-I HKA3A KRTHA3A
Immunogen	Synthesized peptide derived from part region of human protein AA range: 276-326
Specificity	KT33A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Keratin-33B;Keratin;type I cuticular Ha3-II;keratin 33B;K33B;HHA3-II;KT33A;KRT33A HHA3-I HKA3A KRTHA3A;type I cuticular Ha3-I (Hair keratin;type I Ha3-I) (Keratin-33A) (K33A)
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	extracellular space,intermediate filament,
Tissue Specificity	Expressed in the hair follicles.
Function	miscellaneous:There are two types of hair/microfibrillar keratin, I (acidic) and II (neutral to basic).,similarity:Belongs to the intermediate filament family.,tissue specificity:Expressed in the hair follicle.,
Background	This gene encodes a member of the keratin gene family. This gene is one of multiple type I hair keratin genes that are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. As a type I hair keratin, the encoded protein is an acidic protein which heterodimerizes with type II keratins to form hair and nails. There are two isoforms of this protein, encoded by two separate genes, keratin 33A and keratin 33B. [provided by RefSeq, May



2012],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night