



TBD Polyclonal Antibody

Catalog No	YP-Ab-06660
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TUBD1 TUBD
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	TBD Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Delta tubulin;TUBD;TUBD1;Tubulin delta chain;tubulin;delta 1;TBD;TUBD1 TUBD;Tubulin delta chain (Delta-tubulin)
Observed Band	44 kda
Calculated Molecular Weight	51 kDa
Cell Pathway	Nucleus . Cytoplasm . Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Cell projection, cilium . Associated with centrioles. Both cytoplasmic and nuclear. In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head. .
Tissue Specificity	Embryo, Eye,
Function	alternative products:Experimental confirmation may be lacking for some isoforms,function:In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head., similarity:Belongs to the tubulin family.,subcellular location:Associated with centrioles. Both cytoplasmic and nuclear.,

Background

alternative products:Experimental confirmation may be lacking for some isoforms,function:In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head., similarity:Belongs to the tubulin family.,subcellular location:Associated with centrioles. Both cytoplasmic and nuclear.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images