



ZDH13 Polyclonal Antibody

Catalog No	YP-Ab-06693
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ZDHHC13 HIP14L HIP3RP
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZDH13 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DHHC 13;EC:2.3.1.225;HIP14L;HIP14-related protein;HIP3RP;ZDH13; ZDHHC13 HIP14L HIP3RP;Palmitoyltransferase ZDHHC13 (EC 2.3.1.-) (Huntingtin-interacting protein 14-related protein) (HIP14-related protein) (Huntingtin-interacting protein HIP3RP) (Putative MAPK-activating protein PM03) (P
Observed Band	45-50 kDa
Calculated Molecular Weight	622 aa, 71 kDa
Cell Pathway	Golgi apparatus membrane ; Multi-pass membrane protein . Cytoplasmic vesicle membrane ; Multi-pass membrane protein . Low extracellular Mg(2+) induces increase in Golgi and in post-Golgi vesicles. .
Tissue Specificity	Eye,Lung fibroblast,Placenta,Teratocarcinoma,Testis,
Function	catalytic activity:Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA., domain:The DHHC domain is required for palmitoyltransferase activity., function:May be involved in the NF-kappa-B signaling pathway.,sequence caution:Translated as Cys.,similarity:Belongs to the DHHC palmitoyltransferase family. AKR/ZDHHC17 subfamily.,similarity:Contains 1 DHHC-type zinc finger.,similarity:Contains 6 ANK repeats.,

Background

catalytic activity:Palmitoyl-CoA + protein-cysteine = S-palmitoyl protein + CoA.,
domain:The DHHC domain is required for palmitoyltransferase activity.,
function:May be involved in the NF-kappa-B signaling pathway.,sequence
caution:Translated as Cys.,similarity:Belongs to the DHHC palmitoyltransferase
family. AKR/ZDHHC17 subfamily.,similarity:Contains 1 DHHC-type zinc finger.,
similarity:Contains 6 ANK repeats.,

**matters needing
attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For
more information, please consult technical personnel.

Products Images