



B2L13 Polyclonal Antibody

Catalog No	YP-Ab-06722
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	BCL2L13 MIL1 CD003
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	B2L13 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BCL2L13;MIL1;Protein Mil1;B2L13;BCL2L13 MIL1 CD003;Bcl-2-like protein 13 (Bcl2-L-13) (Bcl-rambo) (Protein Mil1)
Observed Band	85 kDa
Calculated Molecular Weight	485 aa, 53 kDa
Cell Pathway	[Isoform 2]: Mitochondrion membrane ; Single-pass membrane protein . Nucleus .; [Isoform 1]: Nucleus .
Tissue Specificity	Ubiquitous, with the highest levels of expression in heart, placenta and pancreas.
Function	alternative products:Experimental confirmation may be lacking for some isoforms,function:May promote the activation of caspase-3 and apoptosis., similarity:Belongs to the Bcl-2 family.,subunit:Monomer.,tissue specificity:Ubiquitous, with the highest levels of expression in heart, placenta and pancreas.,
Background	This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images