



HIPK2 Polyclonal Antibody

Catalog No	YP-Ab-06734
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	HIPK2
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	HIPK2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PRO0593;Homeodomain-interacting protein kinase 2;hHIPk2;FLJ23711; EC:2.7.11.1;HIPK2;Homeodomain-interacting protein kinase 2 (hHIPk2) (EC 2.7.11.1)
Observed Band	131 kDa, 101 kDa
Calculated Molecular Weight	131 kDa
Cell Pathway	Nucleus, PML body . Cytoplasm. Concentrated in PML/POD/ND10 nuclear bodies. Small amounts are cytoplasmic.
Tissue Specificity	Highly expressed in heart, muscle and kidney. Weakly expressed in a ubiquitous way. Down-regulated in several thyroid and breast tumors.
Function	alternative products:Experimental confirmation may be lacking for some isoforms,catalytic activity:ATP + a protein = ADP + a phosphoprotein., function:Protein kinase acting as a corepressor of several transcription factors, including SMAD1 and POU4F1/Brn3a and probably NK homeodomain transcription factors. Inhibits cell growth and promotes apoptosis. Involved in transcriptional activation of TP53 and TP73. Phosphorylation of TP53 may be mediated by a TP53-HIPK2-AXIN1 complex. In response to TGFB, cooperates with DAXX to activate JNK. Phosphorylates the antiapoptotic factor CTBP1 and promotes its proteasomal degradation. In the Wnt/beta-catenin signaling



pathway acts as an intermediate kinase between TAK1 and NLK to promote the proteasomal degradation of MYB (By similarity). Phosphorylates CBX4 upon DNA damage and promotes its E3 SUMO-protein ligase activity.,induction:By UV.,PTM:Phosphory

Background

homeodomain interacting protein kinase 2(HIPK2) Homo sapiens This gene encodes a conserved serine/threonine kinase that is a member of the homeodomain-interacting protein kinase family. The encoded protein interacts with homeodomain transcription factors and many other transcription factors such as p53, and can function as both a corepressor and a coactivator depending on the transcription factor and its subcellular localization. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images