




SG494 Polyclonal Antibody

Catalog No	YP-Ab-06758
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SGK494
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SG494 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ25006;SGK494;SG494;Uncharacterized serine/threonine-protein kinase SgK494 (EC 2.7.11.1) (Sugen kinase 494)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	
Tissue Specificity	Amygdala,Brain,Cerebellum,
Function	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family., similarity:Contains 1 protein kinase domain.,
Background	catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily. Ser/Thr protein kinase family., similarity:Contains 1 protein kinase domain.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night