



RGMB Polyclonal Antibody

Catalog No	YP-Ab-06800
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	RGMB
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	RGMB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RGMB;RGM domain family member B (DRG11-responsive axonal guidance and outgrowth of neurite) (DRAGON)
Observed Band	48kD
Calculated Molecular Weight	48kD
Cell Pathway	Cell membrane ; Lipid-anchor, GPI-anchor . Membrane raft .
Tissue Specificity	Placenta,Testis,
Function	similarity:Belongs to the repulsive guidance molecule (RGM) family.,
Background	RGMB is a glycosylphosphatidylinositol (GPI)-anchored member of the repulsive guidance molecule family (see RGMA, MIM 607362) and contributes to the patterning of the developing nervous system (Samad et al., 2005 [PubMed 15671031]).[supplied by OMIM, Apr 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images