




PAR16 Polyclonal Antibody

Catalog No	YP-Ab-06803
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PARP16 C15orf30
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PAR16 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PAR16;PARP16 C15orf30;Mono [ADP-ribose] polymerase PARP16 (EC 2.4.2.31) (ADP-ribosyltransferase diphtheria toxin-like 15) (ARTD15) (Poly [ADP-ribose] polymerase 16) (PARP-16)
Observed Band	35kD
Calculated Molecular Weight	35kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type IV membrane protein .
Tissue Specificity	Brain,Carcinoma,
Function	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
Background	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC12 cells,
primary antibody was diluted at 1:1000, 4° over night