



SKI Polyclonal Antibody

Catalog No	YP-Ab-06806
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SKI
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SKI Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SKI;Ski oncogene (Proto-oncogene c-Ski)
Observed Band	80kD
Calculated Molecular Weight	80kD
Cell Pathway	Nucleus.
Tissue Specificity	
Function	function:May play a role in terminal differentiation of skeletal muscle cells but not in the determination of cells to the myogenic lineage.,similarity:Belongs to the SKI family.,subunit:Interacts with SMAD2, SMAD3 and SMAD4. Interacts with HIPK2. Part of a complex with HIPK2 and SMAD1/2/3.,
Background	SKI proto-oncogene(SKI) Homo sapiens This gene encodes the nuclear protooncogene protein homolog of avian sarcoma viral (v-ski) oncogene. It functions as a repressor of TGF-beta signaling, and may play a role in neural tube development and muscle differentiation. [provided by RefSeq, Oct 2009],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images