




VMAC Polyclonal Antibody

Catalog No	YP-Ab-06811
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	VMAC
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	VMAC Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VMAC;Vimentin-type intermediate filament-associated coiled-coil protein
Observed Band	18kD
Calculated Molecular Weight	18kD
Cell Pathway	Cytoplasm . Colocalizes with vimentin-type intermediate filaments. .
Tissue Specificity	Cervix,Thymus,
Function	subcellular location:Colocalizes with vimentin-type intermediate filaments.,
Background	subcellular location:Colocalizes with vimentin-type intermediate filaments.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night