



LIN41 Polyclonal Antibody

Catalog No	YP-Ab-06813
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	TRIM71 LIN41
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LIN41 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LIN41;TRIM71 LIN41;Tripartite motif-containing protein 71 (Protein lin-41 homolog)
Observed Band	95kD
Calculated Molecular Weight	95kD
Cell Pathway	Cytoplasm, P-body .
Tissue Specificity	Specifically expressed in testis.
Function	function:May be involved in controlling the timing of organ formation during development.,induction:Down-regulated by retinoic acid in Tera-2 cells.,similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 filamin repeat.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 B box-type zinc fingers.,similarity:Contains 6 NHL repeats.,tissue specificity:Specifically expressed in testis.,
Background	The protein encoded by this gene is an E3 ubiquitin-protein ligase that binds with miRNAs and maintains the growth and upkeep of embryonic stem cells. This gene also is involved in the G1-S phase transition of the cell cycle. [provided by RefSeq, Dec 2015],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images