



LRIG3 Polyclonal Antibody

Catalog No	YP-Ab-06836
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	LRIG3 LIG3 UNQ287/PRO326/PRO335
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LRIG3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA3016;LIG 3;LIG3;LRIG3;LRIG3 LIG3 UNQ287/PRO326/PRO335;Leucine-rich repeats and immunoglobulin-like domains protein 3 (LIG-3)
Observed Band	117 kDa
Calculated Molecular Weight	123 kDa
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Cytoplasmic vesicle membrane ; Single-pass type I membrane protein . Detected in cytoplasmic vesicles when coexpressed with ERBB4. .
Tissue Specificity	Widely expressed.
Function	similarity:Contains 17 LRR (leucine-rich) repeats.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Widely expressed.,
Background	similarity:Contains 17 LRR (leucine-rich) repeats.,similarity:Contains 3 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Widely expressed.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

