



ISLR2 Polyclonal Antibody

Catalog No	YP-Ab-06856
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ISLR2 KIAA1465 LINX UNQ1885/PRO4329
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ISLR2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	241227D6;Immunoglobulin superfamily containing leucine-rich repeat protein 2; KIAA1465;Leucine-rich repeat domain and immunoglobulin domain-containing axon extension protein;LINX;ISLR2;ISLR2 KIAA1465 LINX UNQ1885/PRO4329; Immunoglobulin superfamily containing leucine-rich repeat protein 2 (Leucine-rich repeat domain and immunoglobulin domain-containing axon extension protein)
Observed Band	79 kDa, 130 kDa
Calculated Molecular Weight	745 aa, 79 kDa
Cell Pathway	Cell membrane ; Single-pass membrane protein .
Tissue Specificity	Brain,PCR rescued clones,
Function	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,similarity:Contains 5 LRR (leucine-rich) repeats.,
Background	similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,similarity:Contains 5 LRR (leucine-rich) repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

