



FKBP4 Polyclonal Antibody

Catalog No	YP-Ab-06857
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	FKBP4 FKBP52
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	FKBP4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FKBP4;51 kDa FK506-binding protein;52 kDa FK506-binding protein;EC:5.2.1.8; FK506-binding protein 4;FKBP4 FKBP52;Peptidyl-prolyl cis-trans isomerase FKBP4 (PPIase FKBP4) (EC 5.2.1.8) (51 kDa FK506-binding protein) (FKBP51) (52 kDa FK506-binding protein) (52 kDa FKBP) (FKBP-52) (59 kDa immunophilin) (p59) (FK506-
Observed Band	52 kDa
Calculated Molecular Weight	52 kDa
Cell Pathway	Cytoplasm, cytosol . Mitochondrion . Nucleus . Cytoplasm, cytoskeleton . Cell projection, axon . Shuttles from mitochondria to nucleus; co-localizes in mitochondria with the glucocorticoid receptor (PubMed:21730050). Colocalized with MAPT/TAU in the distal part of the primary cortical neurons (By similarity). .
Tissue Specificity	Widely expressed.
Function	catalytic activity:Peptidylproline (omega=180) = peptidylproline (omega=0)., function:Component of unactivated mammalian steroid receptor complexes that sediment at 8-10 S. May have a rotamase activity. May play a role in the intracellular trafficking of heterooligomeric forms of steroid hormone receptors., PTM:Phosphorylation by CK2 results in loss of HSP90 binding activity (By similarity). Phosphorylated upon DNA damage, probably by ATM or ATR.,



sequence caution:Wrong choice of frame.,similarity:Contains 2 PPlase FKBP-type domains.,similarity:Contains 3 TPR repeats.,subunit:Interacts with NR3C1 and dynein (By similarity). Associates with HSP90 and HSP70 in unactivated steroid hormone receptor complexes. Also interacts with peroxisomal phytanoyl-CoA alpha-hydroxylase (PHYH). Interacts with HSF1 in the HSP90 complex., tissue specificity:Widely expressed.,

Background

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds to the immunosuppressants FK506 and rapamycin. It has high structural and functional similarity to FK506-binding protein 1A (FKBP1A), but unlike FKBP1A, this protein does not have immunosuppressant activity when complexed with FK506. It interacts with interferon regulatory factor-4 and plays an important role in immunoregulatory gene expression in B and T lymphocytes. This encoded protein is known to associate with phytanoyl-CoA alpha-hydroxylase. It can also associate with two heat shock proteins (hsp90 and hsp70) and thus may play a role in the intracellular trafficking of hetero-oligomeric forms of the steroid hormone receptors. This protein corr

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images