



SFRP4 Polyclonal Antibody

Catalog No	YP-Ab-06860
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	SFRP4 FRPHE
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SFRP4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FRP 4;FRPHE;sFRP 4;SFRP4;SFRP4 FRPHE;Secreted frizzled-related protein 4 (sFRP-4) (Frizzled protein;human endometrium) (FrpHE)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Secreted . Cytoplasmic in ovarian tumor cells.
Tissue Specificity	Expressed in mesenchymal cells. Highly expressed in the stroma of proliferative endometrium. Expressed in cardiomyocytes. Shows moderate to strong expression in ovarian tumors with expression increasing as the tumor stage increases. In ovarian tumors, expression levels are inversely correlated with expression of CTNNB1 (at protein level).
Function	disease:Associated with tumor-induced osteomalacia, a disorder in which there is an increase in renal phosphate excretion and a reduction in serum phosphate levels leading to hyperphosphaturia, hypophosphatemia and rickets., domain:The FZ domain is involved in binding with Wnt ligands.,function:Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types. SFRP4 may act as a regulator of adult uterine morphology and function. Increases apoptosis during ovulation possibly



through modulation of FZ1/FZ4/WNT4 signaling (By similarity). Has phosphaturic effects by specifically inhibiting sodium-dependent phosphate uptake.,induction:Increased levels in failing myocardium. Up-regulated in several tumor types including ostomalacia-associated tumors and en

Background

Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night