




OLIG3 Polyclonal Antibody

Catalog No	YP-Ab-06861
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	OLIG3 BHLHB7 BHLHE20
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	OLIG3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OLIG3;OLIG3 BHLHB7 BHLHE20;Oligodendrocyte transcription factor 3 (Oligo3) (Class B basic helix-loop-helix protein 7) (bHLHb7) (Class E basic helix-loop-helix protein 20) (bHLHe20)
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Nucleus .
Tissue Specificity	Ovary,
Function	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
Background	function:May determine the distinct specification program of class A neurons in the dorsal part of the spinal cord and suppress specification of class B neurons., similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night