




FAM3C Polyclonal Antibody

Catalog No	YP-Ab-06869
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	FAM3C ILEI GS3786
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	FAM3C Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM3C;GS3786;ILEI;Interleukin like EMT inducer;Protein FAM3C;FAM3C ILEI GS3786;Protein FAM3C (Interleukin-like EMT inducer)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Secreted . Cytoplasmic vesicle . Cytoplasmic in some cancer cells.
Tissue Specificity	Present in most secretory epithelia (at protein level).
Function	similarity:Belongs to the FAM3 family.,tissue specificity:Ubiquitous.,
Background	This gene is a member of the family with sequence similarity 3 (FAM3) family and encodes a secreted protein with a GG domain. A change in expression of this protein has been noted in pancreatic cancer-derived cells. [provided by RefSeq, Mar 2010],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night