



PARP6 Polyclonal Antibody

Catalog No	YP-Ab-06878
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	PARP6
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PARP6 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PARP6;Poly [ADP-ribose] polymerase 6 (PARP-6) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 17) (ARTD17)
Observed Band	69kD
Calculated Molecular Weight	69kD
Cell Pathway	
Tissue Specificity	Coronary artery,Embryo,Kidney,PCR rescued clones,PNS,Terato
Function	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
Background	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

