




NRARP Polyclonal Antibody

Catalog No	YP-Ab-06881
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	NRARP
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50
Specificity	NRARP Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Notch-regulated ankyrin repeat-containing protein;NRARP
Observed Band	12kD
Calculated Molecular Weight	12kD
Cell Pathway	
Tissue Specificity	Eye,
Function	function:May play a role in the formation of somites.,similarity:Belongs to the NRARP family.,similarity:Contains 2 ANK repeats.,
Background	function:May play a role in the formation of somites.,similarity:Belongs to the NRARP family.,similarity:Contains 2 ANK repeats.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night