



CLUL1 Polyclonal Antibody

Catalog No	YP-Ab-06883
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	CLUL1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CLUL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Retinal-specific clusterin-like protein;Clusterin-like protein 1;Clusterin like protein 1;clusterin like 1 (retinal);CLUL1;Clusterin-like protein 1 (Retinal-specific clusterin-like protein)
Observed Band	70 kDa
Calculated Molecular Weight	54 kDa
Cell Pathway	Secreted .
Tissue Specificity	Brain,Neuroretina,
Function	similarity:Belongs to the clusterin family.,
Background	similarity:Belongs to the clusterin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

