




C1QT1 Polyclonal Antibody

Catalog No	YP-Ab-06886
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	C1QTNF1 CTRP1 UNQ310/PRO353
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	C1QT1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C1QTNF1;Complement C1q tumor necrosis factor-related protein 1;G protein-coupled receptor-interacting protein;GIP;UNQ310/PRO353;C1QT1;C1QTNF1 CTRP1 UNQ310/PRO353;Complement C1q tumor necrosis factor-related protein 1 (G protein-coupled receptor-interacting protein) (GIP)
Observed Band	30kD
Calculated Molecular Weight	30kD
Cell Pathway	Secreted .
Tissue Specificity	Colon,Skin,
Function	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
Background	similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night