



PRDX5 Polyclonal Antibody

Catalog No	YP-Ab-06894
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	PRDX5 ACR1 SBBI10
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PRDX5 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ACR1;Alu corepressor 1;Antioxidant enzyme B166;AOEB166;AOPP;B166;Liver tissue 2D page spot 71B;Peroxiredoxin 5;mitochondrial;Peroxiredoxin V; Peroxisomal antioxidant enzyme;PLP;PMP20;PRDX5;Prx V;PRX5;PRXV; SBBI10;Thioredoxin peroxidase PMP20;Thioredoxin reductase;TPx type VI; PRDX5 ACR1 SBBI10;Peroxiredoxin-5;mitochondrial (EC 1.11.1.15) (Alu corepressor 1) (Antioxidant enzyme B166) (AOEB166) (Liver tissue 2D-page spot 71B) (PLP) (Peroxiredoxin V) (Prx-V) (Peroxisomal antioxidant enzyme) (
Observed Band	23kD
Calculated Molecular Weight	23kD
Cell Pathway	[Isoform Mitochondrial]: Mitochondrion .; [Isoform Cytoplasmic+peroxisomal]: Cytoplasm . Peroxisome matrix . Imported into peroxisomes via peroxisomal targeting signal 1 receptor PEX5. .
Tissue Specificity	Widely expressed.
Function	catalytic activity:2 R'-SH + ROOH = R'-S-S-R' + H(2)O + ROH.,function:Reduces hydrogen peroxide and alkyl hydroperoxides with reducing equivalents provided through the thioredoxin system. Involved in intracellular redox signaling., similarity:Belongs to the peroxiredoxin 2 family.,similarity:Contains 1 thioredoxin domain.,subunit:Monomer.,tissue specificity:Widely expressed.,



Background

This gene encodes a member of the peroxiredoxin family of antioxidant enzymes, which reduce hydrogen peroxide and alkyl hydroperoxides. The encoded protein may play an antioxidant protective role in different tissues under normal conditions and during inflammatory processes. This protein interacts with peroxisome receptor 1. The crystal structure of this protein in its reduced form has been resolved to 1.5 angstrom resolution. This gene uses alternate in-frame translation initiation sites to generate mitochondrial or peroxisomal/cytoplasmic forms. Three transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night