




PAR11 Polyclonal Antibody

Catalog No	YP-Ab-06898
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	PARP11 C12orf6
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	PAR11 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADP-ribosyltransferase diphtheria toxin-like 11;ARTD11;C12orf6;EC:2.4.2.-; MIB006;PAR11;PARP11 C12orf6;Poly [ADP-ribose] polymerase 11 (PARP-11) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 11) (ARTD11)
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Nucleus, nuclear pore complex . Colocalizes with NUP153 at nuclear pores. .
Tissue Specificity	Bone marrow,Brain,Human uterus endothel primary cell cultur
Function	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain., similarity:Contains 1 WWE domain.,
Background	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,similarity:Contains 1 PARP catalytic domain., similarity:Contains 1 WWE domain.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night