



LIGO2 Polyclonal Antibody

Catalog No	YP-Ab-06903
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	LINGO2 LERN3 LRRN6C UNQ9234/PRO31993
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LIGO2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LIGO2;LINGO2 LERN3 LRRN6C UNQ9234/PRO31993;Leucine-rich repeat and immunoglobulin-like domain-containing nogo receptor-interacting protein 2 (Leucine-rich repeat neuronal protein 3) (Leucine-rich repeat neuronal protein 6C)
Observed Band	66kD
Calculated Molecular Weight	66kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Brain,
Function	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 13 LRR (leucine-rich) repeats.,
Background	similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain., similarity:Contains 13 LRR (leucine-rich) repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

