



SPRL1 Polyclonal Antibody

Catalog No	YP-Ab-06912
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	SPARCL1
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SPRL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Hevin;MAST 9;PIG33;SC1;SPARC like 1 (hevin);SPARC like protein 1; SPARCL1;SPRL1;SPARC-like protein 1 (High endothelial venule protein) (Hevin) (MAST 9)
Observed Band	130-140 kDa
Calculated Molecular Weight	664 aa, 75 kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix .
Tissue Specificity	Highly expressed in lymph node, brain, heart, lung, skeletal muscle, ovary, small intestine, and colon, with lower levels in placenta, pancreas, testis, spleen, and thymus, and no expression in kidney, liver, and peripheral blood leukocytes.
Function	similarity:Belongs to the SPARC family.,similarity:Contains 1 EF-hand domain., similarity:Contains 1 follistatin-like domain.,similarity:Contains 1 Kazal-like domain.,tissue specificity:Highly expressed in lymph node, brain, heart, lung, skeletal muscle, ovary, small intestine, and colon, with lower levels in placenta, pancreas, testis, spleen, and thymus, and no expression in kidney, liver, and peripheral blood leukocytes.,

Background

similarity:Belongs to the SPARC family.,similarity:Contains 1 EF-hand domain.,similarity:Contains 1 follistatin-like domain.,similarity:Contains 1 Kazal-like domain.,tissue specificity:Highly expressed in lymph node, brain, heart, lung, skeletal muscle, ovary, small intestine, and colon, with lower levels in placenta, pancreas, testis, spleen, and thymus, and no expression in kidney, liver, and peripheral blood leukocytes.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images