



FKB1B Polyclonal Antibody

Catalog No	YP-Ab-06927
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	FKBP1B FKBP12.6 FKBP1L FKBP9 OTK4
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50
Specificity	FKB1B Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	12.6 kDa FK506 binding protein;12.6 kDa FKBP;FK506 binding protein 1B; FKBP 12.6;FKBP 1B;FKBP12.6;FKBP1B;FKBP1L;FKBP9;h FKBP 12; Immunophilin FKBP12.6;OTK4;PKBP1L;PPlase;PPlase FKBP1B;Rotamase; FKB1B;FKBP1B FKBP12.6 FKBP1L FKBP9 OTK4;Peptidyl-prolyl cis-trans isomerase FKBP1B (PPlase FKBP1B) (EC 5.2.1.8) (12.6 kDa FK506-binding protein) (12.6 kDa FKBP) (FKBP-12.6) (FK506-binding protein 1B) (FKBP-1B) (Immunophilin FKBP12.6) (Rotam
Observed Band	11kD
Calculated Molecular Weight	11kD
Cell Pathway	Cytoplasm . Sarcoplasmic reticulum .
Tissue Specificity	Detected in heart muscle (at protein level). Isoform 1 and isoform 2 are ubiquitous with highest levels in brain and thymus.
Function	catalytic activity:Peptidylproline (omega=180) = peptidylproline (omega=0)., enzyme regulation:Inhibited by both FK506 and rapamycin.,function:Associates with the ryanodine receptor (RYR-2) in cardiac muscle sarcoplasmic reticulum and may play a unique physiological role in excitation-contraction coupling in cardiac muscle. There are four molecules of FKBP12.6 per heart muscle RYR. Has the potential to contribute to the immunosuppressive and toxic effects of



FK506 and rapamycin. PPlases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides., similarity:Belongs to the FKBP-type PPlase family. FKBP1 subfamily., similarity:Contains 1 PPlase FKBP-type domain.,tissue specificity:Ubiquitous for both isoforms with highest levels in brain and thymus.,

Background

The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It is highly similar to the FK506-binding protein 1A. Its physiological role is thought to be in excitation-contraction coupling in cardiac muscle. There are two alternatively spliced transcript variants of this gene encoding different isoforms. [provided by RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night