



LRRN1 Polyclonal Antibody

Catalog No	YP-Ab-06929
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	LRRN1 KIAA1497 Nbla10449 UNQ693/PRO1338
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	LRRN1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FIGLER3;KIAA1497;Leucine-rich repeat neuronal protein 1;Lrrn1; LRRN1_HUMAN;Nbla1449;Neuronal leucine-rich repeat protein 1;NLRR-1; LRRN1 KIAA1497 Nbla10449 UNQ693/PRO1338;Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1)
Observed Band	78kD
Calculated Molecular Weight	73 kDa; 81kD(Calculated).
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Brain,Neuroblastoma,Testis,
Function	similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 11 LRR (leucine-rich) repeats.,
Background	similarity:Contains 1 fibronectin type-III domain.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Contains 11 LRR (leucine-rich) repeats.,
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images