



ZN366 Polyclonal Antibody

Catalog No	YP-Ab-07019
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	ZNF366
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ZN366 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN366;ZNF366;Zinc finger protein 366
Observed Band	81kD
Calculated Molecular Weight	81kD
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in immature and mature dendritic cells (DCs). Not detected in other blood cell types.
Function	function:May be involved in transcriptional regulation.,similarity:Contains 11 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Contains 11 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

