




T4S18 Polyclonal Antibody

Catalog No	YP-Ab-07086
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	TM4SF18
Immunogen	Synthesized peptide derived from human protein . at AA range: 20-100
Specificity	T4S18 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T4S18;TM4SF18;Transmembrane 4 L6 family member 18
Observed Band	22kD
Calculated Molecular Weight	22kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Kidney,
Function	similarity:Belongs to the L6 tetraspanin family.,
Background	similarity:Belongs to the L6 tetraspanin family.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HCT116 cells,
primary antibody was diluted at 1:1000, 4° over night