




FETUB Polyclonal Antibody

Catalog No	YP-Ab-07092
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	FETUB
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Specificity	FETUB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FETUB;16G2;Gugu;Fetuin-like protein IRL685;FetuinB;Fetuin-B (16G2) (Fetuin-like protein IRL685) (Gugu)
Observed Band	42kD
Calculated Molecular Weight	42kD
Cell Pathway	Secreted .
Tissue Specificity	Liver and testis.
Function	similarity:Belongs to the fetuin family.,similarity:Contains 2 cystatin domains., tissue specificity:Liver and testis.,
Background	fetuin B(FETUB) Homo sapiens The protein encoded by this gene is a member of the fetuin family, part of the cystatin superfamily of cysteine protease inhibitors. Fetuins have been implicated in several diverse functions, including osteogenesis and bone resorption, regulation of the insulin and hepatocyte growth factor receptors, and response to systemic inflammation. This protein may be secreted by cells. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night