



LDLR Polyclonal Antibody

Catalog No	YP-Ab-07098
Isotype	IgG
Reactivity	Human;Rat;Mouse;Bovine
Applications	WB;ELISA
Gene Name	LDLR
Immunogen	Synthesized peptide derived from human protein . at AA range: 540-620
Specificity	LDLR Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1A2A3;FH;FHC;LDL receptor;LDLR;Low-density lipoprotein receptor (LDL receptor)
Observed Band	95 kDa
Calculated Molecular Weight	95 kDa
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Membrane, clathrin-coated pit . Golgi apparatus . Early endosome . Late endosome . Lysosome . Rapidly endocytosed upon ligand binding. .
Tissue Specificity	Brain,Liver,Lymph,Plasma,
Function	disease:Defects in LDLR are the cause of familial hypercholesterolemia (FH) [MIM:143890]; a common autosomal semi-dominant disease that affects about 1 in 500 individuals. The receptor defect impairs the catabolism of LDL, and the resultant elevation in plasma LDL-cholesterol promotes deposition of cholesterol in the skin (xanthelasma), tendons (xanthomas), and coronary arteries (atherosclerosis).,function:Binds LDL, the major cholesterol-carrying lipoprotein of plasma, and transports it into cells by endocytosis. In order to be internalized, the receptor-ligand complexes must first cluster into clathrin-coated pits. In case of HIV-1 infection, functions as a receptor for extracellular Tat in neurons, mediating its internalization in uninfected cells.,online information:LDLR

mutation database,online information:The Singapore human mutation and polymorphism database,PTM:N- and O-glycosyla

Background

The low density lipoprotein receptor (LDLR) gene family consists of cell surface proteins involved in receptor-mediated endocytosis of specific ligands. Low density lipoprotein (LDL) is normally bound at the cell membrane and taken into the cell ending up in lysosomes where the protein is degraded and the cholesterol is made available for repression of microsomal enzyme 3-hydroxy-3-methylglutaryl coenzyme A (HMG CoA) reductase, the rate-limiting step in cholesterol synthesis. At the same time, a reciprocal stimulation of cholesterol ester synthesis takes place. Mutations in this gene cause the autosomal dominant disorder, familial hypercholesterolemia. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Sep 2010],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images