




CRIS3 Polyclonal Antibody

Catalog No	YP-Ab-07104
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CRISP3
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	CRIS3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Aeg2;CRISP-3;CRS3;SGP28;dJ442L6.3;CRIS3_HUMAN;CRISP3;Specific granule protein of 28 kDa (SGP28);CRIS3;Cysteine-rich secretory protein 3 (CRISP-3) (Specific granule protein of 28 kDa) (SGP28)
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Secreted. In neutrophils, localized in specific granules.
Tissue Specificity	Salivary gland, pancreas and prostate > epididymis, ovary, thymus and colon.
Function	similarity:Belongs to the CRISP family.,subcellular location:In neutrophils, localized in specific granules.,tissue specificity:Salivary gland, pancreas and prostate > epididymis, ovary, thymus and colon.,
Background	similarity:Belongs to the CRISP family.,subcellular location:In neutrophils, localized in specific granules.,tissue specificity:Salivary gland, pancreas and prostate > epididymis, ovary, thymus and colon.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night