



OTUD4 Polyclonal Antibody

Catalog No	YP-Ab-07111
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	OTUD4 HIN1 KIAA1046
Immunogen	Synthesized peptide derived from human protein . at AA range: 910-990
Specificity	OTUD4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DUBA6;HIN1;HIV 1 induced protein HIN 1;HSHIN1;KIAA1046;OTU;OTUD4; OTUD4 HIN1 KIAA1046;OTU domain-containing protein 4 (HIV-1-induced protein HIN-1)
Observed Band	140-150 kDa
Calculated Molecular Weight	1113 aa, 124 kDa
Cell Pathway	Cytoplasm . Nucleus . Primarily cytoplasmic. .
Tissue Specificity	Brain,Colon carcinoma,Epithelium,Mammary carcinoma,Testis,Tongue,
Function	induction:By HIV-1 insertion.,sequence caution:Contaminating sequence. Originates from a chimeric transcript in HIV-1-infected cells. The premature stop may be due to intron retention.,sequence caution:Intron retention., similarity:Contains 1 OTU domain.,
Background	Alternatively spliced transcript variants have been found for this gene. The smaller protein isoform encoded by the shorter transcript variant is found only in HIV-1 infected cells. [provided by RefSeq, Jul 2010],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images