




ARC Polyclonal Antibody

Catalog No	YP-Ab-07144
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	ARC KIAA0278
Immunogen	Synthesized peptide derived from human protein . at AA range: 70-150
Specificity	ARC Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hArc;Arg3.1;ARC/ARG3.1;Activity-regulated gene 3.1 protein homolog;Activity-regulated cytoskeleton-associated protein;ARC;ARC KIAA0278;Activity-regulated cytoskeleton-associated protein (ARC/ARG3.1) (Activity-regulated gene 3.1 protein homolog) (Arg3.1)
Observed Band	43kD
Calculated Molecular Weight	43kD
Cell Pathway	Extracellular vesicle membrane ; Lipid-anchor . Cell junction, synapse, postsynaptic cell membrane ; Lipid-anchor . Cell junction, synapse . Cell junction, synapse, postsynaptic density . Early endosome membrane . Cell projection, dendrite . Cytoplasm, cytoskeleton . Cytoplasm, cell cortex . Cell projection, dendritic spine . Cytoplasmic vesicle, secretory vesicle, acrosome . Forms virion-like extracellular vesicles that are released from neurons. Enriched in postsynaptic density of dendritic spines. Targeted to inactive synapses following interaction with CAMK2B in the kinase inactive state. Accumulation at weaker synapses may be required to prevent their undesired enhancement. Associated with the cell cortex of neuronal soma and dendrites (By similarity). Associated with the sperm tail (
Tissue Specificity	Brain,Forebrain,Uterus,



Function	function:Required for consolidation of synaptic plasticity as well as formation of long-term memory. Regulates endocytosis of AMPA receptors in response to synaptic activity. Required for homeostatic synaptic scaling of AMPA receptors., similarity:Belongs to the ARC/ARG3.1 family.,subcellular location:Associated with the cell cortex of neuronal soma and dendrites. Enriched in postsynaptic density of dendritic spines. Associated with the sperm tail.,subunit:Interacts with SH3GL1/endophilin-2, SH3GL3/endophilin-3 and DNM2/DYN2.,
Background	function:Required for consolidation of synaptic plasticity as well as formation of long-term memory. Regulates endocytosis of AMPA receptors in response to synaptic activity. Required for homeostatic synaptic scaling of AMPA receptors., similarity:Belongs to the ARC/ARG3.1 family.,subcellular location:Associated with the cell cortex of neuronal soma and dendrites. Enriched in postsynaptic density of dendritic spines. Associated with the sperm tail.,subunit:Interacts with SH3GL1/endophilin-2, SH3GL3/endophilin-3 and DNM2/DYN2.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night