



ZSC10 Polyclonal Antibody

Catalog No	YP-Ab-07163
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ZSCAN10 ZNF206
Immunogen	Synthesized peptide derived from human protein . at AA range: 160-240
Specificity	ZSC10 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	OFNS;ZFP26;ZNF26;Zkscan1;ZSC1_HUMAN;ZSCAN1;Zinc finger protein 26; Q2D61_MOUSE;ZSC10;ZSCAN10 ZNF206;Zinc finger and SCAN domain-containing protein 10 (Zinc finger protein 206)
Observed Band	79kD
Calculated Molecular Weight	8 kDa
Cell Pathway	Nucleus .
Tissue Specificity	PCR rescued clones,
Function	function:May be involved in transcriptional regulation.,similarity:Contains 1 SCAN box domain.,similarity:Contains 14 C2H2-type zinc fingers.,
Background	function:May be involved in transcriptional regulation.,similarity:Contains 1 SCAN box domain.,similarity:Contains 14 C2H2-type zinc fingers.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

