




CDX4 Polyclonal Antibody

Catalog No	YP-Ab-07182
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CDX4
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	CDX4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	caudal type homeobox 4;Caudal type homeobox protein 4;CDX4;Homeobox protein CDX 4;Homeobox protein CDX-4 (Caudal-type homeobox protein 4)
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Belongs to the Caudal homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,
Background	This gene encodes a member of a small subfamily of homeobox containing transcription factors. The encoded protein may regulate homeobox gene expression during anteroposterior patterning and hematopoiesis. [provided by RefSeq, Aug 2012],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night