




MYCL1 Polyclonal Antibody

Catalog No	YP-Ab-07184
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	MYCL1 BHLHE38 LMYC MYCL
Immunogen	Synthesized peptide derived from human protein . at AA range: 170-250
Specificity	MYCL1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein L-Myc;MYCL;LMYC;Class E basic helix-loop-helix protein 38;bHLHe38; MYCL1;MYCL1 BHLHE38 LMYC MYCL;Protein L-Myc-1 (Class E basic helix-loop-helix protein 38) (bHLHe38)
Observed Band	40kD
Calculated Molecular Weight	40kD
Cell Pathway	Nucleus .
Tissue Specificity	Placenta,
Function	similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Binds DNA as an heterodimer with MAX.,
Background	similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein. Binds DNA as an heterodimer with MAX.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night