




CRIS2 Polyclonal Antibody

Catalog No	YP-Ab-07200
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CRISP2 GAPDL5 TPX1
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	CRIS2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRISP-2;CT36;GAPDL5;TPX1;TSP1;CRIS2_HUMAN;CRISP2;Cancer/testis antigen 36 (CT36);Testis-specific protein TPX-1;CRIS2;CRISP2 GAPDL5 TPX1; Cysteine-rich secretory protein 2 (CRISP-2) (Cancer/testis antigen 36) (CT36) (Testis-specific protein TPX-1)
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Secreted .
Tissue Specificity	Testis and epididymis.
Function	function:May regulate some ion channels' activity and thereby regulate calcium fluxes during sperm capacitation.,similarity:Belongs to the CRISP family.,tissue specificity:Testis and epididymis.,
Background	function:May regulate some ion channels' activity and thereby regulate calcium fluxes during sperm capacitation.,similarity:Belongs to the CRISP family.,tissue specificity:Testis and epididymis.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night