




CREG2 Polyclonal Antibody

Catalog No	YP-Ab-07202
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CREG2
Immunogen	Synthesized peptide derived from human protein . at AA range: 30-110
Specificity	CREG2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein CREG2;CREG2
Observed Band	31kD
Calculated Molecular Weight	31kD
Cell Pathway	Secreted .
Tissue Specificity	Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.
Function	PTM:It is not sure whether N-glycosylation is on Asn-165 and/or Asn-166., similarity:Belongs to the CREG family.,tissue specificity:Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.,
Background	PTM:It is not sure whether N-glycosylation is on Asn-165 and/or Asn-166., similarity:Belongs to the CREG family.,tissue specificity:Brain specific mainly in the limbic system and faintly in the spinal cord but not in cerebellum.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night