




C1QT4 Polyclonal Antibody

Catalog No	YP-Ab-07240
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	C1QTNF4 CTRP4
Immunogen	Synthesized peptide derived from human protein . at AA range: 40-120
Specificity	C1QT4 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTRP4;ZACRP4;711E1Rik;9434J15Rik;Adin;C1QT4_HUMAN;C1QTNF4;C1q/TNF-related protein 4;Q4ZJN5_MOUSE;C1QT4_MOUSE;C1QT4;C1QTNF4 CTRP4;Complement C1q tumor necrosis factor-related protein 4
Observed Band	36kD
Calculated Molecular Weight	36kD
Cell Pathway	Secreted .
Tissue Specificity	Widely expressed at low levels (PubMed:21658842). Highest levels in adipocyte tissue and brain (PubMed:24366864).
Function	similarity:Contains 2 C1q domains.,
Background	similarity:Contains 2 C1q domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night